



## Flat Specimen Holder 670.1

The Flat Specimen Holder 670.1 is the standard choice for measuring powder samples under normal laboratory conditions. It is ideal for the generation of very fine peaks at high intensities.

The powder is applied to a  $6\mu m$  thick mylar mounting foil  $10 \times 20mm$  in size. The X-ray assymmetrically penetrates the sample surface at an angle of  $45^\circ$ , whilst the sample oscillates in its plane at approximately 1 Hz with an amplitude of 10mm. This compensates for intensity variations caused by the variously orientated powder particles.

