

## 2-Circle Goniometer 414/414A

The 2-Circle Goniometer 414/414A basically consists of two Goniometers 410 mounted one on top of the other.

The rotational axis of the lower circle (Theta) is fed through the upper circle (2-Theta) by a drive shaft.


This system, when mounted on a detector arm, is ideal for application as an analyser in 4- and 6-Circle Diffractometers.

- four precision configurations
- two load configurations

### Application specific versions:


- vacuum suitable
- antimagnetic
- radiation resistant
- in black

### Precision configurations:



	X1	X2	X3	XE
Accuracy ["]:	<=40	30	15	on request
Repeatability (unidir.) ["]:	<=4	2	2	on request
Reversal error ["]:	<=20	10	7	on request
Eccentricity [ $\mu\text{m}$ ]:	<=5	3	2	2
Wobble ["]:	<=8	4	3	3

### Dimensions [mm]:



D:	L:	H:
138	171	145

## Accessories:

---

Motors:	2-3-/5-Ph.
	Servo/DC
Hand wheels:	0056
Gear boxes:	2056.05
	2056.10
	2056.20
Limit switch 2-Theta-circle:	included
Zero-point control:	9100
Encoder <b>XE</b> :	incremental absolute
Control system:	9300
Goniometer head mounting:	1412, 1413
Arms:	analyser, counter weight
Counter weight:	optional
Base plate:	423.12 included

## Maximum load Theta circle (vertical goniometer axis):

---

Max. load [N] <b>W1</b> :	1450
Max. load [N] <b>W2</b> :	2300

## Maximum load 2-Theta circle:

---

- vertical goniometer axis (V- **W1**, V- **W2**)
- horizontal goniometer axis (H- **W1**, H- **W2**)
- The maximum load on the 2-Theta arm is depending on the distance (momentum).

## Specifications:

---

Travel range [°]:	360
Material (housing/worm gear):	Aluminium/Bronze
Gear ratio:	360:1 / 180:1*
Resolution [°]:	0.001 / 0.002**
	0.00005 / 0.0001* minimum (with gear ratio 20:1)

Min. drive torque [Nm]: 0.10 / 0.13\*

Weight [kg]: 13

Flange size [mm]: 56

\* values concerning 2-Circle Goniometer 414A

\*\* step motor, 1000 steps/revolution

