

## Eulerian Cradle 512.5

### General information:

The Eulerian cradles of the series 500 can be combined with the goniometers of the series 400 to create multi-circle diffractometers. These can be used for analytical investigations in the fields of X-ray and neutron diffraction.

The Eulerian cradle is a full-circle cradle with an asymmetrical design. The Phi- and Chi-circle planes are at right angles to one another.

The motor and signal currents are transmitted via slip rings. This enables an unrestricted rotation around the Phi-axis.

A manual Z-adjustment 5104.A05 is incorporated in the Phi-circle, enabling goniometer heads of the series 1000 to be mounted. Also available is a motorised Z-adjustment.

For sample adjustment an optical microscope or an optional CCD-camera is integrated (see Accessories MiniVID).

Both circles are equipped with zero-point controls and step motors. A range of different motor types and specifications is available according to customer requirements.

### Specifications:

Sphere of confusion [mm]:	0.04***	
Parallelity (Chi-plane to Phi-axis) ["]:	<= +/- 20***	
Weight [kg]:	60	
	Phi-circle	Chi-circle
Travel range [°]:	360	160
Gear ratio:	360:1 / 180:1**	360:1
Accuracy ["]:	30	30
Repeatability (unidir.) ["]:	<= 2	<= 3
Reversal error ["]:	<= 10	<= 15
Resolution [°]:	0.001* / 0.002* **	0.001*
Min. drive torque [Nm]:	0.10 / 0.13**	1.5
Flange size [mm]:	56	83


\* step motor, 1000 steps/revolution

\*\* using Goniometer 410A

\*\*\* with a load of 10kg

## Dimensions [mm]:

---



D:

400

Z:

63.5

## Accessories:

---

Motors:	included
Limit switches:	included
Zero-point control:	included
Gear boxes:	2056.05 (Phi)
	2056.10 (Phi)
	2056.20 (Phi)
	2083.05 (Chi)
	2083.10 (Chi)
Encoder:	2083.20 (Chi)
	incremental
	absolute
Control system:	9300
CCD-camera:	MiniVID
Z-adjustment motorised:	5104.A05M
Goniometer heads:	Series 1000

